

# 1. Uromic JIVE - examination couch with inbuilt urodynamic unit



Modern diagnostic couch used also for interventions in all patient positions: sitting, lying, gynecological use, or standing (male voiding position).

X-ray permeable material that enables its usage for videourodynamics with a RTG C arm. Optimized for urodynamics and suitable also for cystoscopy, ultrasonography, transrectal prostate punctures or urogynecological examinations.

Length of components: seat 400 mm, backrest 500 mm, headrest 400 mm, additional feet 550 mm.

Couch is fully equipped: fixation for accessories (holder for profilometer, cystometric pump and halogen spotlight). Removable additional gynecologic footrests and handrests. Height adjustable electrically between 640 – 940 mm, angle of the backrest adjustable electrically from 0° to +80°, manual adjustment of the backrest angle  $\pm 30^\circ$ , seat (Trendelenburg) adjustable electrically from 0° to +30°. Equipped with electronics; handset control of the 3 motors.

Optional accessories for the examination couch:



a. Urodynamic unit Uromic Jive



b. Ultrasound module EchoJive



c. Uroflow transducer on mobile trolley



d. Cystometry pump hung on the frame



e. Transrectal/transvaginal transducer



f. Profilometry puller on adjustable arm