

No: 10150-000

ERBEJET® 2

A new dimension in waterjet surgery: The ERBEJET 2 in the VIO System.



- Simple operation via the display
- Easy to integrate, compatible with existing systems
- Can be combined with the VIO electrosurgical system
- No thermal trauma to the margins of the incision, the surrounding tissue is left intact
- Extremely uniform water jet; delicate, precise and reproducible effects
- Integrated suction with individually adjustable run-on time
- Ergonomically designed single-use instruments
- Foot switch with ReMode function: remote control straight from the operating table
- Individually configurable programs in simple text: switch on the unit and begin working
- Spares blood vessels, nerves and organs

Technical Data

Mains connection	
System voltage	120 - 240 V
Power frequency	50 / 60 Hz

Dimensions and weight	
Weight	11 kg

Additional product information	
Communication Bus	{ERBE Communication Bus (ECB)}
Packaging unit	1 Piece

Standards	
Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	CF