

No: 10140-100

## VIO® 300 D

The entire electrosurgical spectrum in a modular single unit.



- Automatic output dosage for all regulative technologies
- Precise power output adjusted to the individual indication
- New and improved CUT and COAG functions
- Variable configuration with socket modules
- Upgrades for specific medical specialties are possible at any time
- Simple, interactive and safe operation using the TFT display and FocusView
- Individually configurable programs in simple text: switch on the unit and begin working
- Plug and play: plug in the instrument, start work
- ReMode function: "remote control" via the handle or foot switch straight from the sterile field
- Assignment of the footswitch: the display makes it possible to assign the footswitch or even just one footswitch pedal to the output socket of choice

### Technical Data

<b>Power output</b>	
Maximum CUT output	300 watt at 500 ohm
Safety system	NESSY
Frequency	350 kHz
Maximum COAG output	up to 200 watt

<b>Mains connection</b>	
Power frequency	50 / 60 Hz
Power input during stand-by	40 Watts
Power input during max RF output	500 watt / 920 VA
System voltage	100-120 / 220-240 V
Potential equalization connection	yes

<b>Dimensions and weight</b>	
Width x Height x Depth	410 x 165 x 380 mm
Weight	9,5 kg

<b>Additional product information</b>	
Packaging unit	1 Piece
Display	TFT-Display and FocusView

<b>Standards</b>	
Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	CF
Category acc. to UMDNS	11-490