

No: 10134-000

APC 2: Argon-Plasma Coagulation

Effective hemostasis and devitalization with Argon Plasma Coagulation (APC).



- Non-contact procedure, no sticking of the instrument tip to the tissue
- Effective, even surface coagulation, for uniform hemostasis and devitalization
- Better dosage of penetration depth leads to safer application
- Minimal carbonization and smoke plume formation, minimal pollution through unpleasant odors
- Safe procedure with few complications

Technical Data

Dimensions and weight	
Width x Height x Depth	410 x 80 x 370 mm
Weight	4,8 kg

Gas supply	
Gas type	Argon 4.8 (99.998 %)

Additional product information	
Packaging unit	1 Piece

Standards	
Classification acc. to 93/42 EWG	IIb
Type acc. to EN 60 601-1	CF
Protection category acc. to EN 60 601-1	