

No: 10448-000



ERBOKRYO® CA

Cryosurgery - Versatile and Effective.



- Numerous indications treatable using wide selection of cryoprobes and probe tips, useful in many surgical specialties
- Many interventions can be carried out on an outpatient basis without anesthesia
- No bleeding during treatment or when the necrotized tissue is shed
- Almost no scarring after healing
- Coolant: nitrous oxide (N₂O) or carbon dioxide (CO₂)
- All cryoprobes can be steam autoclaved up to 134° C (short cycle)

Technical Data

Mains connection	
System voltage	230 V
Power frequency	50 / 60 Hz

Dimensions and weight	
Width x Height x Depth	327 x 960 x 400 mm
Weight	13,2 kg

Additional product information	
Packaging unit	1 Piece

Standards	
Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	BF