No: 10140-500



VIO® 100 C

The VIO compact unit for surgical outpatient practices



- Intuitive operation
- Maximum CUT output 100 watt
- Finest regulation possible in the lower power range (1-20 watt) due to digital display, enabling precise operation with microsurgical needles
- Configurable sockets (monopolar, bipolar, neutral)
- AUTO START
- 4 Predefined programs (overwriteable)
- Applications: interventions performed on an outpatient basis in the fields of dermatology, plastic surgery, ENT, gynecology, casualty surgery, ophthalmology

Technical Data

Power output	
Maximum CUT output	up to 100 watt
Maximum COAG output	up to 80 watt
Safety system	yes
Frequency	510-550 kHz

Mains connection	
System voltage	100-120 / 220-240 V
Power frequency	50 / 60 Hz
Power input during stand-by	15 Watts
Power input during max RF output	{190 Watt / 200 VA}
Potential equalization connection	yes

Dimensions and weight	
Width x Height x Depth	280 x 135 x 300 mm
Weight	4 kg

Additional product information			
Packaging unit	1 Piece		

Standards		
Classification acc. to 93/42 EWG	IIb	
Protection category acc. to EN 60 601-1	1	
Type acc. to EN 60 601-1	CF	
Category acc. to UMDNS		