

No: 10142-000

## NT 2: nerve test

Monopolar electrostimulator for the intraoperative localization of nerves and muscle structures



- Versatile application for all surgical procedures in sensitive structures
- Nerves are protected
- The post-operative result is optimized

## Technical Data

### Mains connection

Potential equalization connection	yes
-----------------------------------	-----

### Dimensions and weight

Width x Height x Depth	410 x 80 x 370 mm
Weight	4,05 kg

### Additional product information

Packaging unit	1 Piece
----------------	---------

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	CF