

No: 10322-000

IES 2: Smoke Plume Evacuation

Intelligent smoke plume evacuation with IES 2: Automatically adjustable, quiet and effective for a better environment in the OR; reduced risk of viral and bacterial dispersion.



- High flow rate yet quiet operation
- Activation can be carried out either manually, by foot switch or automatically/synchronized to the electro-surgical application
- Optimized accessories for different applications
- Connection to the VIO System via ERBE Communication Bus (ECB)
- Highly cost-efficient due to long filter life

Technical Data

| Mains connection | |
|-----------------------------------|------------|
| Power frequency | 50 / 60 Hz |
| Potential equalization connection | yes |

| Dimensions and weight | |
|------------------------|--------------------|
| Width x Height x Depth | 205 x 260 x 404 mm |
| Weight | 10,7 kg |

| Additional product information | |
|--------------------------------|---------|
| Packaging unit | 1 Piece |

| Standards | |
|---|----|
| Classification acc. to 93/42 EWG | I |
| Protection category acc. to EN 60 601-1 | I |
| Type acc. to EN 60 601-1 | CF |