

# Technical Specification Vacuson 60LP

---

<b>Model</b>	<b>Vacuson 60LP</b>
Code No.	4179-230 / 4179-115
Input voltage	230V - 50 Hz / 115V - 60Hz
Power consumption	250 VA
Control unit	Vacuum adjustable with rotary knob Manometer for indicating vacuum value On / Off - Switch
Noise level of the pump	Max. 60 dB(A)
Air flow rate	60 l / min
Vacuum (max.)	- 0.9 bar / 686 mmHG
Pump type	piston pump / oil free
Jar capacity	2 liters
Jar	Polysulphone / autoclaveable at max.135°C / transparent, not colored
Jar lid	Polysulphone / autoclaveable at max.135°C / with overflow protection
Vacuum meter	yes ( 0 .... -1.0 bar / class 2.5 )
Suction control valve	yes ( rotary valve )
Dimensions control unit (HxLxD)	260x320x270 mm
Weight of the control unit	13 kg
Approvals	CE 1275, Class 2a IEC 60601, Class1, Type BF

## **Infiltration Pump:**

Flow rate	2 to 22 l / h
Pressure	0 to 2.5 bar
Power consumption	9 Watt

## **Conform Canula Handpiece:**

Stroke	2mm
Frequency	3500 strokes / min.
Motor	E-Motor 31ESS / autoclaveable max. 135°C
Power consumption	100 Watt
Motor cable	4m