

No: 10141-000

## VEM 2: VIO® extension module

The VIO extension module has two additional sockets allowing you an even greater choice of sockets.



- 2 additional sockets can be connected to the VIO electro-surgical unit
- As a result, more flexibility in selecting socket assignment

## Technical Data

### Additional product information

Packaging unit	1 Piece
----------------	---------

### Standards

Classification acc. to 93/42 EWG	IIb
----------------------------------	-----

Protection category acc. to EN 60 601-1	
---	--

Type acc. to EN 60 601-1	CF
--------------------------	----