

No: 20195-310

BiCision® S

The model for open surgery procedures



- Larger thermofusion zone
- Precise dissection exactly in the center
- Homogenous, smooth cut
- Rapid cooling of the jaws
- Optimal preparation characteristics
- Pi benefit: The bowl shaped jaw creates a wider fusion zone than any other instrument with a 5 mm shaft. This makes for additional safety.

Technical Data

Dimensions and weight	
Shaft diameter	5 mm
Length	200 mm

Basic features	
Compatibility	VIO 300 D
Connection	MF plug

Cleaning and Reprocessing	
Sterile packaging	EO
Single-use product	yes

Standards	
Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information	
Packaging unit	5 Piece
Max. electrical capacity	0.25 kVp