

No: 20195-213

BiClamp® 200 C

e.g. for intestinal surgery, stomach surgery or abdominal hysterectomies



- Tissue structures can be coagulated over larger areas
- The individual dissection of vessels is usually not necessary
- Effective and rapid coagulation due to the optimal synchronization of the instrument and the ERBE VIO generator
- Highly economical: ERBE BiClamp is reusable; the procedure saves time and suture material
- AUTO STOP is user-friendly and offers optimal results
- When the instrument is plugged in, Plug & Play automatically adjusts the specific parameters

Technical Data

Dimensions and weight	
Length	200 mm
Cable length	with connecting cable 4 m

Basic features	
Shape / Forming	bent 18°
Connection	MF plug

Surface	
Branches	smooth
Coating	with ceramic coating

Cleaning and Reprocessing	
Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Standards	
Classification acc. to 93/42 EWG	IIB
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information	
Max. electrical capacity	0.25 kVp
Packaging unit	1 Piece