

## ArthVIEW® Telescopes for large joint, 4.0mm, 175mm length



**FULL HD**

● 0° ● 30° ● 45° ● 70°

Colour coded rings identify the „direction of view“

Removable eyepiece for C-mount connection or RZ QuickSnap locking

3-step light post adaptors to connect the most common light cables directly to the endoscope



Sapphire glasses ensure highest stability against damages on the lens system

Optimised glass fibre arrangements ensure an uniform image brightness of highest quality


Full HD high definition pictures with increased detail-resolution through newest CAD-designed rod-lens systems

Fully Autoclavable at 134°C and 2.3bar


Laser welded metal parts




| 225-140-000 ArthVIEW telescope 0° direction of view, 4.0mm diameter, 175mm length



| 225-140-030 Same, 30° direction of view



| 225-140-045 Same, 45° direction of view



| 225-140-070 Same, 70° direction of view