

No: 10402-000

## ERBECRYO® 2

ERBECRYO 2: For Biopsies, Recanalization and Devitalization in Bronchoscopy



### Technical Data

Mains connection	
System voltage	100-120 / 220-240 V
Power frequency	50 / 60 Hz
Potential equalization connection	yes

Dimensions and weight	
Width x Height x Depth	410 X 130 X 370 mm
Weight	6,7 kg

Additional product information	
Packaging unit	1 Piece
Programmspeicherplätze	[10]

Standards	
Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	CF