

## VIO® 3

VIO 3 is yet another of our milestones in technology following VIO 300/200. Take advantage of the progress that you and your patients can now leverage using VIO 3.



- Logical, convenient and visual operator guidance (stepGUIDE)
- Reliable, reproducible tissue effects thanks to state-of-the-art processor technology
- ReMode: up to 6 different settings selectable from the sterile field
- Sockets: easy replacement, connection recommendations, slot assignment, up to 6 instruments freely selectable

## Technical Data

Power output	
Safety system	NESSY

Mains connection	
System voltage	100-120 / 220-240 V
Power frequency	50 / 60 Hz
Power input during stand-by	<30 Watts
Power input during max RF output	[550 Watt]
Potential equalization connection	yes

Dimensions and weight	
Width x Height x Depth	415 x 215 x 375 mm
Weight	12 kg

Additional product information	
Packaging unit	1 Piece

Standards	
Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	CF