

No: 10322-000

IES 2: Smoke Plume Evacuation

Intelligent smoke plume evacuation with IES 2: Automatically adjustable, quiet and effective for a better environment in the OR; reduced risk of viral and bacterial dispersion.



- High flow rate yet quiet operation
- Activation can be carried out either manually, by foot switch or automatically/synchronized to the electrosurgical application
- Optimized accessories for different applications
- Connection to the VIO System via ERBE Communication Bus (ECB)
- Highly cost-efficient due to long filter life

Technical Data

Mains connection	
Power frequency	50 / 60 Hz
Potential equalization connection	yes

Dimensions and weight	
Width x Height x Depth	205 x 260 x 404 mm
Weight	10,7 kg

Additional product information	
Packaging unit	1 Piece

Standards	
Classification acc. to 93/42 EWG	I
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	CF