


RZ Medizintechnik GmbH  
Unter Hasslen 20 | D-78532 Tuttlingen  
Phone: +49 (0)7462/9470-0  
Fax: +49 (0)7462/9470-50  
sales@rz-medizintechnik.com  
www.rz-medizintechnik.com



RZ  
Medizintechnik  
GmbH

*Progress in Surgery*

## CARDIAC Scissors



news  
08.2013

## THE RZ CARDIAC SCISSORS

The new CARDIAC SCISSORS are an improvement of existing and established scissors used in cardiac surgery. In precision and delicacy, the scissors offer a greatly improved and enhanced cutting feel.

These extremely precise and extra fine scissors do have a soft touch and allow finest cutting of vessels with their needle-like scissors tip and a very sharp blades.

In addition, the Super Cut version offers a micro serrated blade preventing the fine vessels slipping away as the scissor cuts.



Different angles available



25°



45°




60°



90°



125°



**Diethrich**  
delicate  
fine long blades  
14mm blade length

**180mm**

125-175-025	25°
125-175-045	45°
125-175-060	60°
125-175-090	90°
125-175-125	125°



**Diethrich**  
extra delicate  
extra fine long blades  
14mm blade length

**180mm**


125-170-025	25°
125-170-045	45°
125-170-060	60°
125-170-090	90°
125-170-125	125°



**Diethrich Hegemann**  
standard blades  
13mm blade length

**180mm**

125-180-025	25°
125-180-045	45°
125-180-060	60°
125-180-090	90°
125-180-125	125°



**SUPER CUT**

126-180-025	25°
126-180-045	45°
126-180-060	60°
126-180-090	90°
126-180-125	125°



**Diethrich Peripheral**  
standard blades  
13mm blade length


**130mm**

125-130-025	25°
125-130-045	45°
125-130-060	60°
125-130-090	90°
125-130-125	125°



**SUPER CUT**

126-130-025	25°
126-130-045	45°
126-130-060	60°
126-130-090	90°
126-130-125	125°



**Diethrich**  
delicate, short blades  
extremely smooth-running  
9mm blade length

**180mm**

125-178-025	25°
125-178-045	45°
125-178-060	60°
125-178-090	90°
125-178-125	125°



**SUPER CUT**

126-178-025	25°
126-178-045	45°
126-178-060	60°
126-178-090	90°
126-178-125	125°



Diethrich Favaloro  
standard blades  
13mm blade length



SUPER CUT

180mm

125-185-125

125°

126-185-125

125°



Diethrich Favaloro  
short blades  
very fine  
9mm blade length



SUPER CUT

170mm

125-168-125


125°

126-168-125

125°



**Diethrich Artery**  
standard blade  
with guide  
13mm blade length



**SUPER CUT**

**180mm**

125-200-125


125°

126-200-125

125°



**DeBakey**  
standard blade  
blunt  
26mm blade length



**DeBakey**  
standard blade  
sharp  
26mm blade length

125-216-025	25°
125-216-045	45°
125-216-060	60°

**180mm**

125-217-025	25°
125-217-045	45°
125-217-060	60°

125-220-025	25°
125-220-045	45°
125-220-060	60°

**225mm**

125-221-025	25°
125-221-045	45°
125-221-060	60°

125-280-025	25°
125-280-045	45°
125-280-060	60°

**280mm**

125-281-025	25°
125-281-045	45°
125-281-060	60°



**Potts-Smith**  
22mm blade  
length



**190mm**

125-490-025	25°
125-490-045	45°
125-490-060	60°



**Potts-Smith**  
22mm blade  
length




**140mm**

125-440-025	25°
125-440-045	45°
125-440-060	60°



**Beall**  
Circumflex



**150mm**

125-300-115	115°
-------------	------



- non-reflecting surface
- more durable
- longer resistance

All models are available  
as CeraCut® version

