

RZ Medizintechnik GmbH
Unter Hasslen 20 | D-78532 Tuttlingen

Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50

sales@rz-medizintechnik.com
www.rz-medizintechnik.com



RZ
Medizintechnik
GmbH

Progress in Surgery

MICS CARDIO VASCULAR INSTRUMENTS




news
03.2014

MICS CARDIO VASCULAR

MICS CARDIO VASCULAR

These minimal invasive instruments are characterized by their innovative and application-oriented design. The very low weight reduces symptoms of fatigue during the surgery. Optimal grip and comfort are ensured by the special surface and shape of the so-called golfball handle. The individual colored handles allow an easier differentiation of the instrument types.



Golfball handle




With catch / Without catch




Axial irrigation connector

MICS CARDIO VASCULAR



RED HANDLE

SCISSORS




	Art.-No.	Working length	Ø	Shape
	295-035-045	250mm	3.5mm	Potts Scissors 45° angled, single action
	295-035-090	250mm	3.5mm	Potts Scissors 90° angled, single action
	295-035-125	250mm	3.5mm	Potts Scissors 125° angled, single action
	295-050-000	250mm	5mm	Metzenbaum Scissors straight, single action
	295-050-010	250mm	5mm	Metzenbaum Scissors straight, double action
	295-050-015	250mm	5mm	Metzenbaum Scissors 15° curved, double action
	295-050-020	250mm	5mm	Jameson Scissors 15° curved, double action
	295-050-030	250mm	5mm	Metzenbaum Scissors 30° curved, double action
	295-050-070	250mm	5mm	Metzenbaum Scissors 70° curved, double action
	295-050-125	250mm	5mm	Potts Scissors 125° angled, single action
	295-070-015	250mm	7mm	Metzenbaum Scissors 15° curved, double action

GOLD HANDLE

NEEDLE HOLDER

with catch




	Art.-No.	Working length	Ø	Shape
	295-035-200	250mm	3.5mm	TC Mini Needle Holder straight, 2x6mm
	295-035-201	250mm	3.5mm	TC Mini Needle Holder curved, 2x6mm
	295-050-200	250mm	5mm	TC Needle Holder straight, 3x10mm
	295-050-201	250mm	5mm	TC Needle Holder curved left, 3x10mm
	295-050-202	250mm	5mm	TC Needle Holder curved right, 3x10mm
	295-050-210	250mm	5mm	TC Ryder Needle Holder straight, 1x8mm
	295-050-212	250mm	5mm	TC Ryder Needle Holder straight, 1x10mm
	295-050-214	250mm	5mm	TC Needle Holder straight, 2x8mm
	295-050-216	250mm	5mm	Needle Holder straight, 3x10mm, diamond dusted
	295-050-217	250mm	5mm	Needle Holder curved left, 3x10mm, diamond dusted
	295-050-219	250mm	5mm	TC Crile Wood Needle Holder curved left, 3x10mm



GREEN HANDLE



GRASPING FORCEPS



	Art.-No. without catch	Art.-No. with catch	Working length	∅	Shape
	295-035-253	295-036-253	250mm	3.5mm	DeBakey Forceps 2x1 teeth, 45° lateral angled, double action
	295-050-250	295-051-250	250mm	5mm	DeBakey Forceps 2x1 teeth, 1.5x10mm, double action
	295-050-251	295-051-251	250mm	5mm	DeBakey Forceps 2x1 teeth, 2.8x10mm, double action
	295-050-252	295-051-252	250mm	5mm	DeBakey Forceps 2x1 teeth, 1.5x11mm, single action
	295-050-253	295-051-253	250mm	5mm	DeBakey Forceps 2x1 teeth, 45° angled, 1.5x14mm, single action
	295-050-254	295-051-254	250mm	5mm	DeBakey Forceps 3x2 teeth, 4.8x10mm, double action
	295-050-255	295-051-255	250mm	5mm	DeBakey Forceps 2x1 teeth, 45° lateral angled, 1.5x14mm, single action
	295-050-260	295-051-260	250mm	5mm	Knot Pusher 45° angled, single action
	295-050-265	295-051-265	250mm	5mm	Resano Grasping Forceps straight, single action
	295-050-270	295-051-270	250mm	5mm	Suture Forceps conical, double action
	295-050-280	295-051-280	250mm	5mm	Gemini Forceps curved, smooth, double action

BLUE HANDLE

SPECIAL INSTRUMENTS



Art.-No.	Working length	Ø	Shape
295-040-300	250mm	4mm	Valve Rongeur oval, 4mm, straight, single action
295-040-301	250mm	4mm	Valve Rongeur oval, 4mm, 30° angled, single action

